

Imaging techniques for low-dose X-ray CT based on hybrid pixels : Compressed-sensing, a solution ?

Yannick Boursier - Groupe imXgam

Centre de Physique des Particules de Marseille



Journée du Traitement du Signal et de l'Image pour le Biomédical
10 septembre 2010



I - Biomedical at CPPM

CPPM : a lab from IN2P3 for particle physics.

Physics experiences : Antares, Atlas, LHCb, D0, and imXgam !

Some imXgam projects :

XPIX :

Hybrid pixels for X-ray : XPAD cameras.

PIXSCAN :

Micro CT-Scanner based on hybrid pixels.

ClearPET/XPAD :


Simultaneous PET/CT imaging based on hybrid pixels.



Outline

- 1- Biomedical imaging based on hybrid pixels
- 2 - Tomography reconstruction
- 3 - Future challenges. Sparsity and Compressed Sensing ?
- 4 - Conclusion

Hybrid pixel detectors for High Energy Physics




Hybrid pixel detectors for High Energy Physics

Inner detector


10^8 hybrid pixels
 $400 \times 50 \mu\text{m}^2$

Hybrid pixels

Hybrid pixel




CCD




Fundamental difference with other detectors (CCDs-like) :
Photon counting mode !
No charge integration !

Hybrid pixels


Hybrid pixel



CCD



Fundamental difference with other detectors (CCDs-like) :
Photon counting mode !
No charge integration !




- No Dark noise
- Energy selection
- Very large dynamic range


Journée du Traitement du Signal et de l'Image pour le Biomédical
 10 septembre 2010

Hybrid pixels

Hybrid pixel




CCD



- very fast data acquisition
- choice of du substrat (Si, CdTE, AsGa)
 - No Dark noise
 - Energy selection
 - Very large dynamic range

Fundamental difference with other detectors (CCDs-like) :
Photon counting mode !
No charge integration !



Photon counting versus charge integration



Journée du Traitement du Signal et de l'Image pour le Biomédical
10 septembre 2010



Photon counting versus charge integration

Photon counting mode

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

- ✓ Very low counting rate

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

- ✓ Very low counting rate
(0,01 ph/pixel/s)

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

- ✓ Very low counting rate
(0,01 ph/pixel/s)



Low statistics

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

- ✓ Very low counting rate
(0,01 ph/pixel/s)
- ✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



Low statistics

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

✓ Very low counting rate
(0,01 ph/pixel/s)

✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



Low statistics

**Enhanced detectability
at low contrast**

Energy threshold of detection

Choice of sensor

Photon counting versus charge integration

Photon counting mode

✓ Very low counting rate
(0,01 ph/pixel/s)



Low statistics

✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



**Enhanced detectability
at low contrast**

Energy threshold of detection

✓ Energy selection of X-rays

Choice of sensor

Photon counting versus charge integration

Photon counting mode

✓ Very low counting rate
(0,01 ph/pixel/s)



Low statistics

✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



Enhanced detectability
at low contrast

Energy threshold of detection

✓ Energy selection of X-rays



Discriminating
diffused photons



Enhanced contrast

Choice of sensor

Photon counting versus charge integration

Photon counting mode

✓ Very low counting rate
(0,01 ph/pixel/s)



Low statistics

✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



**Enhanced detectability
at low contrast**

Energy threshold of detection

✓ Energy selection of X-rays



Discriminating
diffused photons



Enhanced contrast

Choice of sensor

✓ High density sensors

Photon counting versus charge integration

Photon counting mode

✓ Very low counting rate
(0,01 ph/pixel/s)



Low statistics

✓ Large dynamic range
(10^{-2} - 10^6 ph/pixel/s , 80dB)



**Enhanced detectability
at low contrast**

Energy threshold of detection

✓ Energy selection of X-rays



Discriminating
diffused photons



Enhanced contrast

Choice of sensor

✓ High density sensors




Better detection
efficiency




**Enhanced detectability
+
Low statistics**

XPAD3 camera : more than 500,000 pixels of 130 μm



New hybrid pixel camera for X-rays XPAD3/Si


- Photon counting
- Silicon sensor : 500 μm thickness.
- 125 x 75 mm^2 : detector size
- 130 x 130 μm^2 : pixel size
- 560 x 960 pixels
- Fast readout and data transfer : up to 200 frames/s (optical fibre and PCIExpress)



Chips 1x1cm assembled in barrettes,
barrettes assembled in tiles.

Whole-body mouse with spatial resolution of 60 μm

Démonstrateur micro-CT PIXSCAN II




Complete system : 3 blocks

OXFORD Instruments X-ray tube

Target Voltage 10 to 90kv, Target Current up to 2 mA

W target, 13 to 40 μm focal spot size, 80 W, 33 degrees Cone Angle

Démonstrateur micro-CT PIXSCAN II




Journée du Traitement du Signal et de l'Image pour le Biomédical
10 septembre 2010

II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.


Image to reconstruct



II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.


Image to reconstruct



II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.


Image to reconstruct



II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.


Image to reconstruct



II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.

Image to reconstruct




II - Tomography reconstruction

- Basis of tomography : data in 1D + angle, object in 2D.



Image to reconstruct

Sinogram




II - Tomography reconstruction

- Basis of tomography : data in 2D + angle, object in 3D.





Sinogram




II - Tomography reconstruction


- Basis of tomography : data in 2D + angle, object in 3D.



Sinogram



Geometry and physics of the problem





The Beer-Lambert law describes the absorption of different materials along the line of sight.

$$I_j(E) = I^0(E) e^{\left(-\int_{r_j} \mu(l, E) dl\right)}$$

Planar detector

Geometry and physics of the problem




The Beer-Lambert law describes the absorption of different materials along the line of sight.

$$I_j(E) = I^0(E) e^{\left(-\int_{r_j} \mu(l, E) dl\right)}$$

Planar detector

Geometry and physics of the problem



The Beer-Lambert law describes the absorption of different materials along the line of sight.

$$I_j(E) = I^0(E) e^{\left(-\int_{r_j} \mu(l, E) dl\right)}$$

Planar detector

Geometry and physics of the problem

A pixel hybrid will measure $y_{j,T} = \int_T^\infty I_j(E) dE$

Under the assumption of a monochromatic beam, we get at energy E_0

$$y_{j,T} = I^0(E_0) e^{-\int_{r_j} \mu(l, E_0) dl}$$

so that
$$\ln \left(\frac{y_{j,T}}{I^0(E_0)} \right) = - \int_{r_j} \mu(l, E_0) dl$$

In a discrete world where the ray r_j cuts all the voxels i , we finally get

$$x_j = \ln \left(\frac{y_{j,T}}{I^0(E_0)} \right) = - \sum_i \alpha_i \mu_i$$

Geometry and physics of the problem

The problem is expressed as a linear problem :

$$X = A\mu$$

where $X = (x_j)_{j \in \mathbb{N}_m}$, $\in \mathbb{R}^m$ are the measures,

and $\mu = (\mu_i)_{i \in \mathbb{N}_n}$, $\in \mathbb{R}^n$ is the attenuation vector to recover

and $A = (\alpha_{i,j})_{(i,j) \in \mathbb{N}_n \times \mathbb{N}_m}$ is the matrix system.

with $n > m$ in general ...

Realistic model incorporating Poisson noise :

$$X = \mathcal{P}(A\mu)$$

Geometry and physics of the problem

State of the art in tomography :

- Filtered Back-Projection based algorithms.

- «Algebraic» Iterative Methods solving a L_2 problem :

$$\hat{\mu} = \arg \min_{\mu} \|X - A\mu\|_2$$

- «Statistical» Iterative Methods incorporating the Poisson noise :

$$\hat{\mu} = \arg \max_{\mu} f(\mu|X) = \arg \max_{\mu} f(X|\mu) f_p(\mu)$$

where $f(\mu|X)$ is the log-likelihood.

First light XPAD3/PIXSCAN II



Reconstruction performed on a GPU AMD/ATI, Algorithm FDK.
But need of 720 projections and $>1\text{mGy/s}$ at 160 mm



Journée du Traitement du Signal et de l'Image pour
10 septembre 2010

III - Future challenges

1 - Reducing the dose may mean :

- Reducing the statistics
- Reducing the number of projections

→ Need to add a priori information

2 - Go towards color imaging !




- Energy selection ! ..
- Acquire directly color information ?



Reducing the statistics

Tests have been performed on a simulated phantom with the software GATE.

openGATE : active collaboration with more than 1000 referenced users.



Two horizontal slices of the cylindrical simulated phantom, with different density balls.

Simulated phantom

Two distinct levels in the phantom with known-density small balls.

Simulations have been performed with :

- Acquired projections of size 100 x 100 pixels, 360 projections.
- A reconstructed volume of 100 x 100 x 100 pixels.
- Activity (Photons/pixels) ranging from 100 to 25600 ph/pix.

Reducing the statistics

Introduction of a priori information in ML-EM algorithm



$$\hat{\mu} = \arg \max_{\mu} f(\mu|X) = \arg \max_{\mu} f(X|\mu) f_p(\mu)$$

Markov prior based on a median filtering

$$f_p(\mu) = C \exp(-\beta U(\mu))$$


$$f_p(\mu_i) = C \exp\left(-\frac{\beta}{2} \sum_b \frac{(\mu_b - m_b)^2}{m_b}\right)$$

Reducing the statistics





Go towards X-color imaging : concept

Contrast magnification






Go towards X-color imaging : concept

Contrast magnification






Go towards X-color imaging : concept

Contrast magnification



Go towards X-color imaging : concept

Capteur CdTe efficace jusqu'à 100 keV



courtesy: A. Butler, MARS Biomedical Imaging Ltd.

3 pharmaceutiques, fantôme perspex



Université Albert-Ludwigs à Freiburg
et université de Canterbury

Go towards X-color imaging : concept

Iode : circulation pulmonaire

Barium : poumons

Calcium : os normaux



Système quantitatif developpé à l'université d'Erlangen-Nürnberg

courtesy: A. Butler, MARS Biomedical Imaging Ltd.

Conclusion

New XPAD3 hybrid pixels camera for X-ray photon counting developed at CPPM :

CT-Scanner based on hybrid pixels.

Simultaneous PET/CT scanner for bimodality images.

Simulations as well as real acquisitions have proven the reliability of hybrid pixels for micro-CT scanner and PET/CT !

Adapted algorithms :

For low dose : Poisson noise to take into account

For small number of projections : regularization needed ... sparsity basis ?

For color imaging : adapt the acquisition framework ... CS theory for help ?

Implementation on GPUs strongly speeds-up the reconstruction (between 100 and 300 times faster compared to CPUs). To be implemented for iterative methods ?

